

EASY

RELIABLE

DURABLE



SAMPLE CHANGER

Positions: 15 beakers of 50 ml

Specially designed for wine analysis

Volume to start titration: 20 ml

Speed: 1 to 2 minutes per sample

Rinse: possible between titrations

Calibration: automated in 2 points

Dimensions (W x D x H) : 53 x 38 x 40 cm

Weight: 5 kg

AUTOMATIC TITRATOR

Supplied ready to use for the determination of titratable acidity and pH

- 25 ml burette
- pH electrode
- Tubing set
- Titration head
- 3 way valve
- Paddle stirrer
- Printer
- 1 liter bottle with drying tube

Easy to use interface

Ability to define and store new programs.

Work end point and inflection point.

Dimensions (W x D x H) : 30 x 25 x 47 cm

Weight: 4 kg